

IL20 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP50996-2**

Capture Antibody Information

Catalog Number: 60691-3-PBS	Clone ID: 1G2G4	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: interleukin 20
Isotype: IgG2b	GenBank: BC074949	Gene ID: 50604
Purification Method: Protein A Magarose purification	Immunogen Catalog Number: Ag12447	

Detection Antibody Information

Catalog Number: 60691-2-PBS	Clone ID: 1F2B12	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: interleukin 20
Isotype: IgG2b	GenBank: BC074949	Gene ID: 50604
Purification Method: Protein A Magarose purification	Immunogen Catalog Number: Ag12447	

Applications

Tested Applications: Cytometric bead array	Range: 1.563-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP50996-2 targets IL20 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: IL20 Monoclonal antibody, PBS Only (Capture) 60691-3-PBS (1G2G4). 100 µg. Concentration 1 mg/mL.

Detection antibody: IL20 Monoclonal antibody, PBS Only (Detector) 60691-2-PBS (1F2B12). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP50996-2, IL20 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60691-3-PBS. Detection antibody: 60691-2-PBS. Standard:Ag12447. Range: 1.563-100 ng/mL.