

CLUAP1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP51001-1

Capture Antibody Information

Catalog Number: 68353-2-PBS	Clone ID: 1A9G4	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: clusterin associated protein 1
Isotype: IgG1	GenBank: BC017070	Gene ID: 23059
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag12722	

Detection Antibody Information

Catalog Number: 68353-3-PBS	Clone ID: 2C11F7	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: clusterin associated protein 1
Isotype: IgG1	GenBank: BC017070	Gene ID: 23059
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag12722	

Applications

Tested Applications: Cytometric bead array	Range: 0.391-50 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP51001-1 targets CLUAP1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CLUAP1 Monoclonal antibody, PBS Only (Capture) 68353-2-PBS (1A9G4). 100 µg. Concentration 1 mg/ml.

Detection antibody: CLUAP1 Monoclonal antibody, PBS Only (Detector) 68353-3-PBS (2C11F7). 100 µg. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP51001-1, CLUAP1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68353-2-PBS. Detection antibody: 68353-3-PBS. Standard: Ag12722. Range: 0.391-50 ng/mL