

For Research Use Only

TYK2 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51012-1

Capture Antibody Information

| | | |
|--|--|--|
| Catalog Number: 67411-2-PBS | Clone ID: 1H5B4 | Conjugate: Unconjugated |
| Host: Mouse | Reactivity: human | Full name: tyrosine kinase 2 |
| Isotype: IgG1 | GenBank: BC014243 | Gene ID: 7297 |
| Purification Method: Protein G Magarose purification | Immunogen Catalog Number: Ag9442 | |

Detection Antibody Information

| | | |
|--|--|--|
| Catalog Number: 67411-3-PBS | Clone ID: 1C10F4 | Conjugate: Unconjugated |
| Host: Mouse | Reactivity: human | Full name: tyrosine kinase 2 |
| Isotype: IgG2a | GenBank: BC014243 | Gene ID: 7297 |
| Purification Method: Protein A Magarose purification | Immunogen Catalog Number: Ag9442 | |

Applications

| | | |
|--|--|---|
| Tested Applications: Cytometric bead array | Range: 0.098-100 ng/mL (Cytometric Bead Array) | Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results. |
|--|--|---|

Product Information

MP51012-1 targets TYK2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TYK2 Monoclonal antibody, PBS Only (Capture) 67411-2-PBS (1H5B4). 100 µg. Concentration 1 mg/ml.

Detection antibody: TYK2 Monoclonal antibody, PBS Only (Detector) 67411-3-PBS (1C10F4). 100 µg. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51012-1, TYK2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67411-2-PBS. Detection antibody: 67411-3-PBS. Standard: Ag9442. Range: 0.098-100 ng/mL.