

For Research Use Only

ERO1L Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51017-2

Capture Antibody Information

Catalog Number: 67416-2-PBS	Clone ID: 1G10C6	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: ERO 1-like (<i>S. cerevisiae</i>)
Isotype: IgG1	GenBank: BC008674	Gene ID: 30001
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag29910	

Detection Antibody Information

Catalog Number: 67416-4-PBS	Clone ID: 2B5E5	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: ERO 1-like (<i>S. cerevisiae</i>)
Isotype: IgG1	GenBank: BC008674	Gene ID: 30001
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag29910	

Applications

Tested Applications: Cytometric bead array	Range: 0.195-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP51017-2 targets ERO1L in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ERO 1L Monoclonal antibody, PBS Only (Capture) 67416-2-PBS (1G10C6). 100 μ g. Concentration 1 mg/ml.

Detection antibody: ERO 1L Monoclonal antibody, PBS Only (Capture/Detector) 67416-4-PBS (2B5E5). 100 μ g. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51017-2, ERO1L Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67416-2-PBS. Detection antibody: 67416-4-PBS. Standard:Ag29910. Range: 0.195-100 ng/mL.