

For Research Use Only

NELL1 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51021-1

Capture Antibody Information

Catalog Number: 60707-1-PBS	Clone ID: 2A9H3	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: NEL-like 1 (chicken)
Isotype: IgG1	GenBank: BC069674	Gene ID: 4745
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag27773	

Detection Antibody Information

Catalog Number: 60707-2-PBS	Clone ID: 2C10G5	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: NEL-like 1 (chicken)
Isotype: IgG1	GenBank: BC069674	Gene ID: 4745
Purification Method: Protein G purification	Immunogen Catalog Number: Ag27773	

Applications

Tested Applications: Cytometric bead array	Range: 0.781-200 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP51021-1 targets NELL1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: NELL1 Monoclonal antibody, PBS Only (Capture) 60707-1-PBS (2A9H3). 100 µg. Concentration 1 mg/ml.

Detection antibody: NELL1 Monoclonal antibody, PBS Only (Detector) 60707-2-PBS (2C10G5). 100 µg. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51021-1, NELL1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60707-1-PBS. Detection antibody: 60707-2-PBS. Standard:Ag27773. Range: 0.781-200 ng/mL.