

For Research Use Only

# RPGRIP1 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51024-2

## Capture Antibody Information

Catalog Number:  
60709-3-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein G purification

Clone ID:  
1B8D6

Reactivity:  
human

GenBank:  
BC039089

Immunogen Catalog Number:  
Ag31294

Conjugate:  
Unconjugated

Full name:  
retinitis pigmentosa GTPase regulator  
interacting protein 1

Gene ID:  
57096

## Detection Antibody Information

Catalog Number:  
60709-2-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein G Magarose purification

Clone ID:  
2G4H10

Reactivity:  
human

GenBank:  
BC039089

Immunogen Catalog Number:  
Ag31294

Conjugate:  
Unconjugated

Full name:  
retinitis pigmentosa GTPase regulator  
interacting protein 1

Gene ID:  
57096

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.391-100 ng/mL (Cytometric Bead  
Array)

Recommended Dilutions:  
It is recommended that this reagent  
should be titrated in each testing  
system to obtain optimal results.

## Product Information

MP51024-2 targets RPGRIP1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: RPGRIP1 Monoclonal antibody, PBS Only (Capture) 60709-3-PBS (1B8D6). 100  $\mu$ g. Concentration 1 mg/ml.

Detection antibody: RPGRIP1 Monoclonal antibody, PBS Only (Detector) 60709-2-PBS (2G4H10). 100  $\mu$ g. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP51024-2, RGRIP1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60709-3-PBS. Detection antibody: 60709-2-PBS. Standard:Ag31294. Range: 0.391-100 ng/mL.