

For Research Use Only

ZP3 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51025-3

Capture Antibody Information

Catalog Number:
60710-1-PBS

Host:
Mouse

Isotype:
IgG2a

Purification Method:
Protein A Magarose purification

Clone ID:
2C10G5

Reactivity:
human

GenBank:
BC113949

Immunogen Catalog Number:
Ag17685

Conjugate:
Unconjugated

Full name:
zona pellucida glycoprotein 3 (sperm receptor)

Gene ID:
7784

Detection Antibody Information

Catalog Number:
60710-3-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
2B10A10

Reactivity:
human

GenBank:
BC113949

Immunogen Catalog Number:
Ag17685

Conjugate:
Unconjugated

Full name:
zona pellucida glycoprotein 3 (sperm receptor)

Gene ID:
7784

Applications

Tested Applications:
Cytometric bead array

Range:
0.391-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51025-3 targets ZP3 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ZP3 Monoclonal antibody, PBS Only (Capture) 60710-1-PBS (2C10G5). 100 μ g. Concentration 1 mg/mL.

Detection antibody: ZP3 Monoclonal antibody, PBS Only (Detector) 60710-3-PBS (2B10A10). 100 μ g. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51025-3, ZP3 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60710-1-PBS. Detection antibody: 60710-3-PBS. Standard:Ag17685. Range: 0.391-100 ng/mL.