

For Research Use Only

JUB Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51044-3

Capture Antibody Information

Catalog Number:

60728-1-PBS

Host:

Mouse

Isotype:

IgG1

Purification Method:

Protein G purification

Clone ID:

2G5C2

Reactivity:

human

GenBank:

BC007580

Immunogen Catalog Number:

Ag29742

Conjugate:

Unconjugated

Full name:

jub, ajuba homolog (Xenopus laevis)

Gene ID:

84962

Detection Antibody Information

Catalog Number:

60728-4-PBS

Host:

Mouse

Isotype:

IgG1

Purification Method:

Protein G Magarose purification

Clone ID:

1B9E7

Reactivity:

human

GenBank:

BC007580

Immunogen Catalog Number:

Ag29742

Conjugate:

Unconjugated

Full name:

jub, ajuba homolog (Xenopus laevis)

Gene ID:

84962

Applications

Tested Applications:

Cytometric bead array

Range:

0.098-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:

It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51044-3 targets JUB in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: JUB Monoclonal antibody, PBS Only (Capture) 60728-1-PBS (2G5C2). 100 µg. Concentration 1 mg/ml.

Detection antibody: JUB Monoclonal antibody, PBS Only (Detector) 60728-4-PBS (1B9E7). 100 µg. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51044-3, JUB Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60728-1-PBS. Detection antibody: 60728-4-PBS. Standard: Ag29742. Range: 0.098-100 ng/mL.