

For Research Use Only

RABL5 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51060-3

Capture Antibody Information

Catalog Number:
60734-5-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
3D8A6
Reactivity:
human
GenBank:
BC038668
Immunogen Catalog Number:
Ag26078

Conjugate:
Unconjugated
Full name:
RAB, member RAS oncogene family-like 5
Gene ID:
64792

Detection Antibody Information

Catalog Number:
60734-6-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
2F3D4
Reactivity:
human
GenBank:
BC038668
Immunogen Catalog Number:
Ag26078

Conjugate:
Unconjugated
Full name:
RAB, member RAS oncogene family-like 5
Gene ID:
64792

Applications

Tested Applications:
Cytometric bead array

Range:
0.098-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51060-3 targets RABL5 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: RABL5 Monoclonal antibody, PBS Only (Capture) 60734-5-PBS (3D8A6). 100 µg. Concentration 1 mg/mL.

Detection antibody: RABL5 Monoclonal antibody, PBS Only (Detector) 60734-6-PBS (2F3D4). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51060-3, RABL5 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60734-5-PBS. Detection antibody: 60734-6-PBS. Standard:Ag26078. Range: 0.098-100 ng/mL.