

For Research Use Only

FUT4 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51063-1

Capture Antibody Information

Catalog Number:
60736-1-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1H3E4

Reactivity:
human

GenBank:
BC136374

Immunogen Catalog Number:
Ag17695

Conjugate:
Unconjugated

Full name:
fucosyltransferase 4 (alpha (1,3) fucosyltransferase, myeloid-specific)

Gene ID:
2526

Detection Antibody Information

Catalog Number:
60736-2-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G purification

Clone ID:
2H9H1

Reactivity:
human

GenBank:
BC136374

Immunogen Catalog Number:
Ag17695

Conjugate:
Unconjugated

Full name:
fucosyltransferase 4 (alpha (1,3) fucosyltransferase, myeloid-specific)

Gene ID:
2526

Applications

Tested Applications:
Cytometric bead array

Range:
0.098-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51063-1 targets FUT4 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: FUT4 Monoclonal antibody, PBS Only (Capture) 60736-1-PBS (1H3E4). 100 μ g. Concentration 1 mg/ml.

Detection antibody: FUT4 Monoclonal antibody, PBS Only (Detector) 60736-2-PBS (2H9H1). 100 μ g. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51063-1, FUT4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60736-1-PBS. Detection antibody: 60736-2-PBS. Standard:Ag17695. Range: 0.098-100 ng/mL.