

For Research Use Only

LRI G1 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51067-3

Capture Antibody Information

Catalog Number:
60738-3-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
3G5A10

Reactivity:
human

GenBank:
BC071561

Immunogen Catalog Number:
Ag17880

Conjugate:
Unconjugated

Full name:
leucine-rich repeats and immunoglobulin-like domains 1

Gene ID:
26018

Detection Antibody Information

Catalog Number:
60738-5-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
2H8E7

Reactivity:
human

GenBank:
BC071561

Immunogen Catalog Number:
Ag17880

Conjugate:
Unconjugated

Full name:
leucine-rich repeats and immunoglobulin-like domains 1

Gene ID:
26018

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51067-3 targets LRI G1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: LRI G1 Monoclonal antibody, PBS Only (Capture) 60738-3-PBS (3G5A10). 100 μ g. Concentration 1 mg/mL.

Detection antibody: LRI G1 Monoclonal antibody, PBS Only (Detector) 60738-5-PBS (2H8E7). 100 μ g. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51067-3, LRIg1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60738-3-PBS. Detection antibody: 60738-5-PBS. Standard:Ag17880. Range: 0.781-100 ng/mL.