

For Research Use Only

WDR78 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51072-3

Capture Antibody Information

Catalog Number: 60742-4-PBS	Clone ID: 2C4A10	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: WD repeat domain 78
Isotype: IgG2a	GenBank: BC126379	Gene ID: 79819
Purification Method: Protein A purification	Immunogen Catalog Number: Ag27018	

Detection Antibody Information

Catalog Number: 60742-5-PBS	Clone ID: 2G2C5	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: WD repeat domain 78
Isotype: IgG1	GenBank: BC126379	Gene ID: 79819
Purification Method: Protein G purification	Immunogen Catalog Number: Ag27018	

Applications

Tested Applications: Cytometric bead array	Range: 0.195-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
------------------------------------------------------	----------------------------------------------------------	-------------------------------------------------------------------------------------------------------------------------------------------

Product Information

MP51072-3 targets WDR78 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: WDR78 Monoclonal antibody, PBS Only (Capture) 60742-4-PBS (2C4A10). 100 µg. Concentration 1 mg/ml.

Detection antibody: WDR78 Monoclonal antibody, PBS Only (Detector) 60742-5-PBS (2G2C5). 100 µg. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51072-3, WDR78 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60742-4-PBS. Detection antibody: 60742-5-PBS. Standard:Ag27018. Range: 0.195-100 ng/mL.