

For Research Use Only

IFT52 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51076-3

Capture Antibody Information

Catalog Number:
60743-1-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1A1A9

Reactivity:
human

GenBank:
BC039831

Immunogen Catalog Number:
Ag11613

Conjugate:
Unconjugated

Full name:
intraflagellar transport 52 homolog
(Chlamydomonas)

Gene ID:
51098

Detection Antibody Information

Catalog Number:
60743-4-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
3D3E6

Reactivity:
human

GenBank:
BC039831

Immunogen Catalog Number:
Ag11613

Conjugate:
Unconjugated

Full name:
intraflagellar transport 52 homolog
(Chlamydomonas)

Gene ID:
51098

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51076-3 targets IFT52 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: IFT52 Monoclonal antibody, PBS Only (Capture) 60743-1-PBS (1A1A9). 100 μ g. Concentration 1 mg/ml.

Detection antibody: IFT52 Monoclonal antibody, PBS Only (Detector) 60743-4-PBS (3D3E6). 100 μ g. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51076-3, IFT52 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60743-1-PBS. Detection antibody: 60743-4-PBS. Standard:Ag11613. Range: 0.781-100 ng/mL.