

For Research Use Only

Cytokeratin 8 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51080-1

Capture Antibody Information

Catalog Number:
68627-1-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G purification

Clone ID:
1F2H5
Reactivity:
human, mouse, rat
GenBank:
BC000654
Immunogen Catalog Number:
Ag0488

Conjugate:
Unconjugated
Full name:
keratin 8
Gene ID:
3856

Detection Antibody Information

Catalog Number:
68627-2-PBS
Host:
Mouse
Isotype:
IgG1
Purification Method:
Protein G Magarose purification

Clone ID:
2A6A3
Reactivity:
human
GenBank:
BC000654
Immunogen Catalog Number:
Ag0488

Conjugate:
Unconjugated
Full name:
keratin 8
Gene ID:
3856

Applications

Tested Applications:
Cytometric bead array

Range:
0.098-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51080-1 targets Cytokeratin 8 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Cytokeratin 8 Monoclonal antibody, PBS Only (Capture) 68627-1-PBS (1F2H5). 100 µg. Concentration 1 mg/mL.

Detection antibody: Cytokeratin 8 Monoclonal antibody, PBS Only (Detector) 68627-2-PBS (2A6A3). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51080-1, Cytokeratin 8 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 68627-1-PBS. Detection antibody: 68627-2-PBS. Standard:Ag0488. Range: 0.098-100 ng/mL