

For Research Use Only

SMARCA4 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51084-1

Capture Antibody Information

Catalog Number:
66561-2-PBS

Host:
Mouse

Isotype:
IgG2b

Purification Method:
Protein A Magarose purification

Clone ID:
2E7B8

Reactivity:
human

GenBank:
BC150298

Immunogen Catalog Number:
Ag16256

Conjugate:
Unconjugated

Full name:
SWI/SNF related, matrix associated, actin dependent regulator of chromatin, subfamily a, member 4

Gene ID:
6597

Detection Antibody Information

Catalog Number:
66561-3-PBS

Host:
Mouse

Isotype:
IgG2a

Purification Method:
Protein A purification

Clone ID:
1F4B10

Reactivity:
human

GenBank:
BC150298

Immunogen Catalog Number:
Ag16256

Conjugate:
Unconjugated

Full name:
SWI/SNF related, matrix associated, actin dependent regulator of chromatin, subfamily a, member 4

Gene ID:
6597

Applications

Tested Applications:
Cytometric bead array

Range:
0.098-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51084-1 targets SMARCA4 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SMARCA4 Monoclonal antibody, PBS Only (Capture) 66561-2-PBS (2E7B8). 100 μ g. Concentration 1 mg/ml.

Detection antibody: SMARCA4 Monoclonal antibody, PBS Only (Detector) 66561-3-PBS (1F4B10). 100 μ g. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51084-1, SMARCA4 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66561-2-PBS. Detection antibody: 66561-3-PBS. Standard: Ag16256. Range: 0.098-100 ng/mL.