

CD268/BAFF-R Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP51091-2

Capture Antibody Information

Catalog Number:
60763-1-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1D6F7

Reactivity:
human

GenBank:
BC105123

Immunogen Catalog Number:
Ag18512

Conjugate:
Unconjugated

Full name:
tumor necrosis factor receptor superfamily, member 13C

Gene ID:
115650

Detection Antibody Information

Catalog Number:
60763-3-PBS

Host:
Mouse

Isotype:
IgM

Purification Method:
Caprylic acid/ammonium sulfate precipitation

Clone ID:
2D9D11

Reactivity:
human

GenBank:
BC105123

Immunogen Catalog Number:
Ag18512

Conjugate:
Unconjugated

Full name:
tumor necrosis factor receptor superfamily, member 13C

Gene ID:
115650

Applications

Tested Applications:
Cytometric bead array

Range:
0.098-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51091-2 targets CD268/BAFF-R in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CD268/BAFF-R Monoclonal antibody, PBS Only (Capture) 60763-1-PBS (1D6F7). 100 µg. Concentration 1 mg/mL.

Detection antibody: CD268/BAFF-R Monoclonal antibody, PBS Only (Detector) 60763-3-PBS (2D9D11). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP51091-2, TNFRSF13C Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60763-1-PBS. Detection antibody: 60763-3-PBS. Standard:Ag18512. Range: 0.098-100 ng/mL