

For Research Use Only

# ATF6 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51101-1

## Capture Antibody Information

Catalog Number:  
66563-2-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein G purification

Clone ID:  
3C3D10

Reactivity:  
human

GenBank:  
BC014969

Immunogen Catalog Number:  
Ag21456

Conjugate:  
Unconjugated

Full name:  
activating transcription factor 6

Gene ID:  
22926

## Detection Antibody Information

Catalog Number:  
66563-3-PBS

Host:  
Mouse

Isotype:  
IgG2b

Purification Method:  
Protein A Magarose purification

Clone ID:  
3A9B12

Reactivity:  
human

GenBank:  
BC014969

Immunogen Catalog Number:  
Ag21456

Conjugate:  
Unconjugated

Full name:  
activating transcription factor 6

Gene ID:  
22926

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.098-50 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP51101-1 targets ATF6 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: ATF6 Monoclonal antibody, PBS Only (Capture) 66563-2-PBS (3C3D10). 100  $\mu$ g. Concentration 1 mg/mL.

Detection antibody: ATF6 Monoclonal antibody, PBS Only (Detector) 66563-3-PBS (3A9B12). 100  $\mu$ g. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage buffer:  
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP51101-1, ATF6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66563-2-PBS. Detection antibody: 66563-3-PBS. Standard:ATF6. Range: 0.098-50 ng/mL