

For Research Use Only

# SMARCA2-Specific Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51103-3

## Capture Antibody Information

Catalog Number:  
60774-4-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein G purification

Clone ID:  
2A2E1

Reactivity:  
human

GenBank:  
NM\_003070

Immunogen Catalog Number:  
Ag33074

Conjugate:  
Unconjugated

Full name:  
SWI/SNF related, matrix associated, actin dependent regulator of chromatin, subfamily a, member 2

Gene ID:  
6595

## Detection Antibody Information

Catalog Number:  
60774-5-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein G Magarose purification

Clone ID:  
1B3A1

Reactivity:  
human

GenBank:  
NM\_003070

Immunogen Catalog Number:  
Ag33074

Conjugate:  
Unconjugated

Full name:  
SWI/SNF related, matrix associated, actin dependent regulator of chromatin, subfamily a, member 2

Gene ID:  
6595

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.391-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP51103-3 targets SMARCA2-Specific in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: SMARCA2-Specific Monoclonal antibody, PBS Only (Capture) 60774-4-PBS (2A2E1). 100 µg. Concentration 1 mg/ml.

Detection antibody: SMARCA2-Specific Monoclonal antibody, PBS Only (Detector) 60774-5-PBS (1B3A1). 100 µg. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:  
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:  
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP51103-3, SMARCA2-Specific Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60774-4-PBS. Detection antibody: 60774-5-PBS. Standard:Ag33074. Range: 0.391-100 ng/mL