

For Research Use Only

# NFKB1 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51119-1

## Capture Antibody Information

Catalog Number:  
66992-2-PBS

Host:  
Mouse

Isotype:  
IgG2b

Purification Method:  
Protein A purification

Clone ID:  
1A1F1

Reactivity:  
human

GenBank:  
BC051765

Immunogen Catalog Number:  
Ag5832

Conjugate:  
Unconjugated

Full name:  
nuclear factor of kappa light polypeptide gene enhancer in B-cells 1

Gene ID:  
4790

## Detection Antibody Information

Catalog Number:  
66992-1-PBS

Host:  
Mouse

Isotype:  
IgG2a

Purification Method:  
Protein A purification

Clone ID:  
2G1E3

Reactivity:  
human, mouse

GenBank:  
BC051765

Immunogen Catalog Number:  
Ag5832

Conjugate:  
Unconjugated

Full name:  
nuclear factor of kappa light polypeptide gene enhancer in B-cells 1

Gene ID:  
4790

## Applications

Tested Applications:  
Cytometric bead array

Range:  
3.125-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP51119-1 targets NFKB1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: NFKB1 Monoclonal antibody, PBS Only (Capture) 66992-2-PBS (1A1F1). 100 µg. Concentration 1 mg/ml.

Detection antibody: NFKB1,p105,p50 Monoclonal antibody, PBS Only (Detector) 66992-1-PBS (2G1E3). 100 µg. Concentration 1 mg/ml.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP51119-1, NFKB1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66992-2-PBS. Detection antibody: 66992-1-PBS. Standard: Ag5832. Range: 3.125-100 ng/mL.