

For Research Use Only

CDH6 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51125-1

Capture Antibody Information

Catalog Number:
60787-1-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
3A4D11

Reactivity:
human

GenBank:
BC000019

Immunogen Catalog Number:
Ag37731

Conjugate:
Unconjugated

Full name:
cadherin 6, type 2, K-cadherin (fetal kidney)

Gene ID:
1004

Detection Antibody Information

Catalog Number:
60787-2-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
2F3D6

Reactivity:
human

GenBank:
BC000019

Immunogen Catalog Number:
Ag37731

Conjugate:
Unconjugated

Full name:
cadherin 6, type 2, K-cadherin (fetal kidney)

Gene ID:
1004

Applications

Tested Applications:
Cytometric bead array

Range:
1.563-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51125-1 targets CDH6 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CDH6 Monoclonal antibody, PBS Only (Capture) 60787-1-PBS (3A4D11). 100 μ g. Concentration 1 mg/mL.

Detection antibody: CDH6 Monoclonal antibody, PBS Only (Detector) 60787-2-PBS (2F3D6). 100 μ g. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51125-1, CDH6 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60787-1-PBS. Detection antibody: 60787-2-PBS. Standard: Ag37731. Range: 1.563-100 ng/mL.