

For Research Use Only

EPB41 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51150-1

Capture Antibody Information

Catalog Number:
60801-1-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1C7C5

Reactivity:
human

GenBank:
BC039079

Immunogen Catalog Number:
Ag4365

Conjugate:
Unconjugated

Full name:
erythrocyte membrane protein band 4.1 (elliptocytosis 1, RH-linked)

Gene ID:
2035

Detection Antibody Information

Catalog Number:
60801-2-PBS

Host:
Mouse

Isotype:
IgG2a

Purification Method:
Protein A purification

Clone ID:
1A12C9

Reactivity:
human

GenBank:
BC039079

Immunogen Catalog Number:
Ag4365

Conjugate:
Unconjugated

Full name:
erythrocyte membrane protein band 4.1 (elliptocytosis 1, RH-linked)

Gene ID:
2035

Applications

Tested Applications:
Cytometric bead array

Range:
3.125-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51150-1 targets EPB41 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: EPB41 Monoclonal antibody, PBS Only (Capture) 60801-1-PBS (1C7C5). 100 μ g. Concentration 1 mg/mL.

Detection antibody: EPB41 Monoclonal antibody, PBS Only (Detector) 60801-2-PBS (1A12C9). 100 μ g. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51150-1, EPB41 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60801-1-PBS. Detection antibody: 60801-2-PBS. Standard:Ag4365. Range: 3.125-100 ng/mL.