

For Research Use Only

p130Cas / BCAR1 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51156-3

Capture Antibody Information

Catalog Number: 67215-5-PBS	Clone ID: 1F3E2	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: breast cancer anti-estrogen resistance 1
Isotype: IgG1	GenBank: BC062556	Gene ID: 9564
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag10604	

Detection Antibody Information

Catalog Number: 67215-6-PBS	Clone ID: 1C12C10	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: breast cancer anti-estrogen resistance 1
Isotype: IgG1	GenBank: BC062556	Gene ID: 9564
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag10604	

Applications

Tested Applications: Cytometric bead array	Range: 6.25-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP51156-3 targets p130Cas / BCAR1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: p130Cas / BCAR1 Monoclonal antibody, PBS Only (Capture) 67215-5-PBS (1F3E2). 100 µg. Concentration 1 mg/mL.

Detection antibody: p130Cas / BCAR1 Monoclonal antibody, PBS Only (Detector) 67215-6-PBS (1C12C10). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51156-3, p130Cas / BCAR1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67215-5-PBS. Detection antibody: 67215-6-PBS. Standard: Ag10604. Range: 6.25-100 ng/mL