

For Research Use Only

CD21 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51164-2

Capture Antibody Information

Catalog Number:
66701-4-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
2D2B11

Reactivity:
human

GenBank:
BC136394

Immunogen Catalog Number:
Ag18058

Conjugate:
Unconjugated

Full name:
complement component (3d/Epstein Barr virus) receptor 2

Gene ID:
1380

Detection Antibody Information

Catalog Number:
66701-5-PBS

Host:
Mouse

Isotype:
IgG2a

Purification Method:
Protein A Magarose purification

Clone ID:
2D10A2

Reactivity:
human

GenBank:
BC136394

Immunogen Catalog Number:
Ag18058

Conjugate:
Unconjugated

Full name:
complement component (3d/Epstein Barr virus) receptor 2

Gene ID:
1380

Applications

Tested Applications:
Cytometric bead array

Range:
0.781-50 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51164-2 targets CD21 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CD21 Monoclonal antibody, PBS Only (Capture) 66701-4-PBS (2D2B11). 100 μ g. Concentration 1 mg/mL.

Detection antibody: CD21 Monoclonal antibody, PBS Only (Detector) 66701-5-PBS (2D10A2). 100 μ g. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51164-2, CD21 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66701-4-PBS. Detection antibody: 66701-5-PBS. Standard: Ag18058. Range: 0.781-50 ng/mL