

PHF11 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number: **MP51172-2**

Capture Antibody Information

Catalog Number: 60808-1-PBS	Clone ID: 1D4D5	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: PHD finger protein 11
Isotype: IgG2a	GenBank: BC017212	Gene ID: 51131
Purification Method: Protein A purification	Immunogen Catalog Number: Ag18041	

Detection Antibody Information

Catalog Number: 60808-3-PBS	Clone ID: 1C12H4	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: PHD finger protein 11
Isotype: IgG2a	GenBank: BC017212	Gene ID: 51131
Purification Method: Protein A Magarose purification	Immunogen Catalog Number: Ag18041	

Applications

Tested Applications: Cytometric bead array	Range: 1.563-200 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
---	---	--

Product Information

MP51172-2 targets PHF11 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: PHF11 Monoclonal antibody, PBS Only (Capture) 60808-1-PBS (1D4D5). 100 µg. Concentration 1 mg/mL.

Detection antibody: PHF11 Monoclonal antibody, PBS Only (Detector) 60808-3-PBS (1C12H4). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP51172-2, PHF11 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60808-1-PBS. Detection antibody: 60808-3-PBS. Standard:Ag18041. Range: 1.563-200 ng/mL.