

For Research Use Only

CCL22/MDC Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51202-2

Capture Antibody Information

Catalog Number: 60754-3-PBS	Clone ID: 2B8G7	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: chemokine (C-C motif) ligand 22
Isotype: IgG2b	GenBank: BC027952	Gene ID: 6367
Purification Method: Protein A Magarose purification	Immunogen Catalog Number: Ag4275	

Detection Antibody Information

Catalog Number: 60754-2-PBS	Clone ID: 3E1F3	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: chemokine (C-C motif) ligand 22
Isotype: IgG1	GenBank: BC027952	Gene ID: 6367
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag4275	

Applications

Tested Applications: Cytometric bead array	Range: 0.781-50 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	---	---

Product Information

MP51202-2 targets CCL22/MDC in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CCL22/MDC Monoclonal antibody, PBS Only (Capture/Detector) 60754-3-PBS (2B8G7). 100 µg. Concentration 1 mg/mL.

Detection antibody: CCL22/MDC Monoclonal antibody, PBS Only (Detector) 60754-2-PBS (3E1F3). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage


Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51202-2, CCL22/MDC Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60754-3-PBS. Detection antibody: 60754-2-PBS. Standard: Ag4275. Range: 0.781-50 ng/mL.