

For Research Use Only

# CCL22/MDC Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51202-4

## Capture Antibody Information

Catalog Number:  
60754-5-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein G purification

Clone ID:  
2G11B7

Reactivity:  
human

GenBank:  
BC027952

Immunogen Catalog Number:  
Ag4275

Conjugate:  
Unconjugated

Full name:  
chemokine (C-C motif) ligand 22

Gene ID:  
6367

## Detection Antibody Information

Catalog Number:  
60754-3-PBS

Host:  
Mouse

Isotype:  
IgG2b

Purification Method:  
Protein A Magarose purification

Clone ID:  
2B8G7

Reactivity:  
human

GenBank:  
BC027952

Immunogen Catalog Number:  
Ag4275

Conjugate:  
Unconjugated

Full name:  
chemokine (C-C motif) ligand 22

Gene ID:  
6367

## Applications

Tested Applications:  
Cytometric bead array

Range:  
0.781-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP51202-4 targets CCL22/MDC in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: CCL22/MDC Monoclonal antibody, PBS Only (Capture) 60754-5-PBS (2G11B7). 100  $\mu$ g. Concentration 1 mg/mL.

Detection antibody: CCL22/MDC Monoclonal antibody, PBS Only (Capture/Detector) 60754-3-PBS (2B8G7). 100  $\mu$ g. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP51202-4, CCL22/MDC Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60754-5-PBS. Detection antibody: 60754-3-PBS. Standard:Ag4275. Range: 0.781-100 ng/mL