

TRIM59 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP51219-2

Capture Antibody Information

Catalog Number: 60824-3-PBS	Clone ID: 1F3C3	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: tripartite motif-containing 59
Isotype: IgG1	GenBank: NM_173084	Gene ID: 286827
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag29572	

Detection Antibody Information

Catalog Number: 60824-4-PBS	Clone ID: 1B5A2	Conjugate: Unconjugated
Host: Mouse	Reactivity: human	Full name: tripartite motif-containing 59
Isotype: IgG1	GenBank: NM_173084	Gene ID: 286827
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag29572	

Applications

Tested Applications: Cytometric bead array	Range: 0.781-100 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
--	--	---

Product Information

MP51219-2 targets TRIM59 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TRIM59 Monoclonal antibody, PBS Only (Capture) 60824-3-PBS (1F3C3). 100 µg. Concentration 1 mg/mL.

Detection antibody: TRIM59 Monoclonal antibody, PBS Only (Detector) 60824-4-PBS (1B5A2). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.


Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

Selected Validation Data



Cytometric bead array standard curve of MP51219-2, TRIM59 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60824-3-PBS. Detection antibody: 60824-4-PBS. Standard:Ag29572. Range: 0.781-100 ng/mL.