

Arfaptin-1 Monoclonal Matched Antibody Pair, PBS Only

Catalog Number:MP51228-1

Capture Antibody Information

Catalog Number:
67461-2-PBS

Host:
Mouse

Isotype:
IgG2a

Purification Method:
Protein A Magarose purification

Clone ID:
1F6G10

Reactivity:
human

GenBank:
BC103759

Immunogen Catalog Number:
Ag12097

Conjugate:
Unconjugated

Full name:
ADP-ribosylation factor interacting protein 1

Gene ID:
27236

Detection Antibody Information

Catalog Number:
67461-3-PBS

Host:
Mouse

Isotype:
IgG2a

Purification Method:
Protein A Magarose purification

Clone ID:
2E3H3

Reactivity:
human

GenBank:
BC103759

Immunogen Catalog Number:
Ag12097

Conjugate:
Unconjugated

Full name:
ADP-ribosylation factor interacting protein 1

Gene ID:
27236

Applications

Tested Applications:
Cytometric bead array

Range:
0.098-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51228-1 targets Arfaptin-1 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: Arfaptin-1 Monoclonal antibody, PBS Only (Capture) 67461-2-PBS (1F6G10). 100 µg. Concentration 1 mg/mL.

Detection antibody: Arfaptin-1 Monoclonal antibody, PBS Only (Capture/Detector) 67461-3-PBS (2E3H3). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:


Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:

PBS only

Selected Validation Data



Cytometric bead array standard curve of MP51228-1, Arfaptin-1 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 67461-2-PBS. Detection antibody: 67461-3-PBS. Standard: Ag12097. Range: 0.098-100 ng/mL.