

For Research Use Only

# TNRC6A Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51246-3

## Capture Antibody Information

Catalog Number:  
60841-4-PBS

Host:  
Mouse

Isotype:  
IgG2a

Purification Method:  
Protein A Magarose purification

Clone ID:  
1B12G12

Reactivity:  
human

GenBank:  
BC024324

Immunogen Catalog Number:  
Ag34199

Conjugate:  
Unconjugated

Full name:  
trinucleotide repeat containing 6A

Gene ID:  
27327

## Detection Antibody Information

Catalog Number:  
60841-5-PBS

Host:  
Mouse

Isotype:  
IgG1

Purification Method:  
Protein G Magarose purification

Clone ID:  
1G3H6

Reactivity:  
human

GenBank:  
BC024324

Immunogen Catalog Number:  
Ag34199

Conjugate:  
Unconjugated

Full name:  
trinucleotide repeat containing 6A

Gene ID:  
27327

## Applications

Tested Applications:  
Cytometric bead array

Range:  
1.563-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:  
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

## Product Information

MP51246-3 targets TNRC6A in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: TNRC6A Monoclonal antibody, PBS Only (Capture) 60841-4-PBS (1B12G12). 100 µg. Concentration 1 mg/mL.

Detection antibody: TNRC6A Monoclonal antibody, PBS Only (Detector) 60841-5-PBS (1G3H6). 100 µg. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

## Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage buffer:

PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP51246-3, TNRC6A Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 60841-4-PBS. Detection antibody: 60841-5-PBS. Standard:Ag34199. Range: 1.563-100 ng/mL.