

For Research Use Only

GM130;GOLGA2 Monoclonal Matched Antibody Pair, PBS Only



Catalog Number:MP51313-3

Capture Antibody Information

Catalog Number:
66662-4-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1G3B2

Reactivity:
human

GenBank:
BC014188

Immunogen Catalog Number:
Ag1848

Conjugate:
Unconjugated

Full name:
golgi autoantigen, golgin subfamily a, 2

Gene ID:
2801

Detection Antibody Information

Catalog Number:
66662-3-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1E12A4

Reactivity:
human

GenBank:
BC014188

Immunogen Catalog Number:
Ag1848

Conjugate:
Unconjugated

Full name:
golgi autoantigen, golgin subfamily a, 2

Gene ID:
2801

Applications

Tested Applications:
Cytometric bead array

Range:
0.391-100 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP51313-3 targets GM130;GOLGA2 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: GM130;GOLGA2 Monoclonal antibody, PBS Only (Capture) 66662-4-PBS (1G3B2). 100 μ g. Concentration 1 mg/mL.

Detection antibody: GM130;GOLGA2 Monoclonal antibody, PBS Only (Detector) 66662-3-PBS (1E12A4). 100 μ g. Concentration 1 mg/mL.

Unconjugated mouse monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP51313-3, GM130;GOLGA2 Monoclonal Matched Antibody Pair, PBS Only. Capture antibody: 66662-4-PBS. Detection antibody: 66662-3-PBS. Standard:Ag1848. Range: 0.391-100 ng/mL