

For Research Use Only

BCL6 Hybrid Matched Antibody Pair, PBS Only



Catalog Number:MP80008-2

Capture Antibody Information

Catalog Number: 66340-3-PBS	Clone ID: 2H4D2	Conjugate: Unconjugated
Host: Mouse	Reactivity: Human	Full name: B-cell CLL/lymphoma 6
Isotype: IgG1	GenBank: BC150184	Gene ID: 604
Purification Method: Protein G Magarose purification	Immunogen Catalog Number: Ag15519	

Detection Antibody Information

Catalog Number: 81756-1-PBS	Clone ID: 2F16	Conjugate: Unconjugated
Host: Rabbit	Reactivity: Human	Full name: B-cell CLL/lymphoma 6
Isotype: IgG	GenBank: BC150184	Gene ID: 604
Purification Method: Protein A purification	Immunogen Catalog Number: Ag15519	

Applications

Tested Applications: Cytometric bead array	Range: 0.156-20 ng/mL (Cytometric Bead Array)	Recommended Dilutions: It is recommended that this reagent should be titrated in each testing system to obtain optimal results.
------------------------------------------------------	---------------------------------------------------------	-------------------------------------------------------------------------------------------------------------------------------------------

Product Information

MP80008-2 targets BCL6 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: BCL6 Monoclonal antibody, PBS Only (Capture) 66340-3-PBS (2H4D2). 100 µg. Concentration 1 mg/ml.

Detection antibody: BCL6 Recombinant antibody, PBS Only (Detector) 81756-1-PBS (2F16). 100 µg. Concentration 1 mg/ml.

Alternative BCL6 matched antibody pairs: MP50215-1, MP80008-1, MP80008-3

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.


Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP80008-2, BCL6 Hybrid Matched Antibody Pair, PBS Only. Capture antibody: 66340-3-PBS. Detection antibody: 81756-1-PBS. Standard: Ag15519. Range: 0.156-20 ng/mL