

For Research Use Only

BCL6 Hybrid Matched Antibody Pair, PBS Only



Catalog Number:MP80008-3

Capture Antibody Information

Catalog Number:
81750-1-PBS

Host:
Rabbit

Isotype:
IgG

Purification Method:
Protein A purification

Clone ID:
6C10

Reactivity:
Human

GenBank:
BC150184

Immunogen Catalog Number:
Ag15519

Conjugate:
Unconjugated

Full name:
B-cell CLL/lymphoma 6

Gene ID:
604

Detection Antibody Information

Catalog Number:
66340-4-PBS

Host:
Mouse

Isotype:
IgG1

Purification Method:
Protein G Magarose purification

Clone ID:
1E5E11

Reactivity:
Human

GenBank:
BC150184

Immunogen Catalog Number:
Ag15519

Conjugate:
Unconjugated

Full name:
B-cell CLL/lymphoma 6

Gene ID:
604

Applications

Tested Applications:
Cytometric bead array

Range:
0.156-20 ng/mL (Cytometric Bead Array)

Recommended Dilutions:
It is recommended that this reagent should be titrated in each testing system to obtain optimal results.

Product Information

MP80008-3 targets BCL6 in immunoassays as a matched antibody pair. Validated in Cytometric bead array.

Capture antibody: BCL6 Recombinant antibody, PBS Only (Capture) 81750-1-PBS (6C10). 100 μ g. Concentration 1 mg/mL.

Detection antibody: BCL6 Monoclonal antibody, PBS Only (Detector) 66340-4-PBS (1E5E11). 100 μ g. Concentration 1 mg/mL.

Alternative BCL6 matched antibody pairs: MP50215-1, MP80008-1, MP80008-2

Unconjugated rabbit recombinant monoclonal antibody pair in PBS only storage buffer at a concentration of 1 mg/mL, ready for conjugation. Created using Proteintech's proprietary in-house recombinant technology. Recombinant production enables unrivalled batch-to-batch consistency, easy scale-up, and future security of supply.

Matched antibody pairs are designed for use in a variety of assays and platforms that require matched antibody pairs.

Antibody use should be optimized for each application and assay.

Storage

Storage:
Store at -80°C.
Storage buffer:
PBS only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP80008-3, BCL6 Hybrid Matched Antibody Pair, PBS Only. Capture antibody: 81750-1-PBS. Detection antibody: 66340-4-PBS. Standard: Ag15519. Range: 0.156-20 ng/mL