

For Research Use Only

# PE Anti-Human CD120a (H398)

Catalog Number: PE-65170



## Basic Information

### Catalog Number:

PE-65170

### Size:

100 tests, 10 µl/test

### Source:

Mouse

### Isotype:

IgG2a

### GenBank Accession Number:

BC010140

### GeneID (NCBI):

7132

### UNIPROT ID:

P19438

### Full Name:

tumor necrosis factor receptor  
superfamily, member 1A

### Calculated MW:

455 aa, 50 kDa

### Purification Method:

The purified antibody is conjugated with R-phycoerythrin (PE) under optimum conditions. The conjugate is purified by size-exclusion chromatography.

### CloneNo.:

H398

### Excitation/Emission maxima wavelengths:

496 nm, 565 nm / 578 nm

## Applications

### Tested Applications:

FC

### Species Specificity:

Human

## Background Information

## Storage

### Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

### Storage Buffer:

PBS with 0.1% sodium azide.

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



100 ul human peripheral blood were surface stained with 10.00 ul PE Anti-Human CD120a (PE-65170, Clone: H398). Cells were not fixed.



100 ul human peripheral blood were surface stained with PE-mouse IgG2a isotype control antibody. Cells were not fixed.