

For Research Use Only

# FcZero-rAb™ PE Anti-Human CD89 Rabbit Recombinant Antibody

Catalog Number: PE-FcA98105



## Basic Information

Catalog Number:

PE-FcA98105

Size:

100 tests, 5 ul/test

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC027953

GeneID (NCBI):

2204

Full Name:

Fc fragment of IgA, receptor for

Calculated MW:

287 aa, 32 kDa

Purification Method:

Protein A purification

CloneNo.:

241426E3

Excitation/Emission maxima  
wavelengths:

496 nm, 565 nm / 578 nm

## Applications

Tested Applications:

FC

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA.

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



1x10<sup>6</sup> human peripheral blood leukocytes were surface stained with 5 ul PE Anti-Human CD89 Rabbit Recombinant Antibody (PE-FcA98105, Clone: 241426E3), or 5 ul PE Rabbit IgG Isotype Control Recombinant Antibody (PE-FcA98136, Clone: 240953C9). Cells were not fixed.