

For Research Use Only.
Not For Use In Diagnostics.

Multi-rAb® Biotin-Goat Anti-Rabbit Recombinant Secondary Antibody (H+L)



Catalog Number: **RGAR009**

Information

Catalog Number:

RGAR009

Host:

Goat

Applications:

WB, IHC, ELISA

Reactivity:

Rabbit

Physical State:

Liquid

Conjugation:

Biotin

Recommended Dilutions

1:100-1:1000 for enzyme immunohistochemistry on tissue sections;
1:5000-1:20000 for ELISA and Western blotting with chromogenic substrates.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 50% glycerol, 10 mg/mL BSA, 0.1% ProClin 300, pH 7.4.

Aliquoting is unnecessary

Purity

The antibody was purified from culture media supernatant by immunoaffinity chromatography using Protein G beads.


For technical support and original validation data for this product please
contact:

T: 4006900926

E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.

Selected Validation Data



Immunohistochemical analysis of paraffin-embedded human appendicitis tissue slide using 80008-1-RR (Alpha smooth muscle actin specific antibody) at dilution of 1:10000. Multi-rAb® Biotin-Goat Anti-Rabbit Recombinant Secondary Antibody (H+L) was used at 4 ug/mL for signal amplification. HRP conjugated streptavidin was used for detection.