

For Research Use Only.
Not For Use In Diagnostics.

HRP-conjugated Affinipure Donkey Anti-Human IgG(H+L)

Catalog Number:SA00001-11



www.ptglab.com

Information

Catalog Number:

SA00001-11

Host:

Donkey

Applications:

ELISA, WB, Dot blot

Reactivity:

Human (Minimal Cross Reactivity:

Bovine, Chicken, Goat, Guinea Pig,

Syrian Hamster, Horse, Rabbit, Mouse,

Rat, Sheep Serum Proteins)

Physical State:

Liquid

Conjugation:

HRP

Recommended Dilutions

1:2000-1:10,000 for western blotting with ECL substrates

1:1000-1:20,000 for ELISA and Western blotting with chromogenic substrates

Stabilizer

3 mg/ml BSA

Safety Notes

This product is for research use only, not for diagnostic or therapeutic use.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

0.01 M Sodium phosphate, 0.25 M NaCl, 50% glycerol, 3 mg/ml BSA, pH 7.6.

Aliquoting is unnecessary for -20°C storage

Purity

The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data