

For Research Use Only.  
Not For Use In Diagnostics.

# HRP-conjugated Affinipure Goat Anti-Guinea pig IgG(H+L)

Catalog Number:SA00001-12



[www.ptglab.com](http://www.ptglab.com)

---

## Information

<b>Catalog Number:</b> SA00001-12	<b>Reactivity:</b> Guinea pig
<b>Host:</b> Goat	<b>Physical State:</b> Liquid
<b>Applications:</b> ELISA, WB, Dot blot	<b>Conjugation:</b> HRP

---

## Recommended Dilutions

1:2000-1:10,000 for western blotting with ECL substrates  
1:1000-1:20,000 for ELISA and Western blotting with chromogenic substrates

---

## Stabilizer

3 mg/ml BSA

---

## Safety Notes

This product is for research use only, not for diagnostic or therapeutic use.

---

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.  
**Storage Buffer:**  
0.01 M Sodium phosphate, 0.25 M NaCl, 50% glycerol, 3 mg/ml BSA, pH 7.6.  
Aliquoting is unnecessary for -20°C storage

---

## Purity

The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.

---

For technical support and original validation data for this product please contact:

T: 4006900926      E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)      W: [ptgcn.com](http://ptgcn.com)

---

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

---

## Selected Validation Data