

For Research Use Only.
Not For Use In Diagnostics.

HRP-conjugated Affinipure Rabbit Anti-Sheep IgG(H+L)

Catalog Number:SA00001-16



www.ptglab.com

Information

Catalog Number:	SA00001-16	Reactivity:	Sheep
Host:	Rabbit	Physical State:	Liquid
Applications:	ELISA, WB, Dot blot	Conjugation:	HRP

Recommended Dilutions

1:2000-1:10,000 for western blotting with ECL substrates
1:1000-1:20,000 for ELISA and Western blotting with chromogenic substrates

Stabilizer

3 mg/ml BSA

Safety Notes

This product is for research use only, not for diagnostic or therapeutic use.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
0.01 M Sodium phosphate, 0.25 M NaCl, 50% glycerol, 3 mg/ml BSA, pH 7.6.
Aliquoting is unnecessary for -20°C storage

Purity

The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.

For technical support and original validation data for this product please contact:

T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data