

For Research Use Only.  
Not For Use In Diagnostics.

# Alkaline Phosphatase–conjugated Affinipure Donkey Anti-Goat IgG(H+L)



[www.ptglab.com](http://www.ptglab.com)

Catalog Number: SA00002-3

## Information

Catalog Number:	SA00002-3	Reactivity:	Goat
Host:	Donkey	Physical State:	Liquid
Applications:	ELISA, WB	Conjugation:	Alkaline phosphatase

## Recommended Dilutions

1:1000 – 1:3000 for ELISA and Western blotting

## Preservative

0.01% Sodium Azide

## Stabilizer

5 mg/ml BSA

## Safety Notes

This product is for research use only, not for diagnostic or therapeutic use.

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.  
**Storage Buffer:**  
0.01 M Tris-HCl, 0.25 M NaCl, 50% glycerol, 5 mg/ml BSA, 0.01% Sodium Azide, pH 8.0.  
Aliquoting is unnecessary for -20°C storage

## Purity

The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.

For technical support and original validation data for this product please contact:

T: 4006900926      E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)      W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

---

## Selected Validation Data