

For Research Use Only.
Not For Use In Diagnostics.

Alkaline Phosphatase–conjugated Affinipure Donkey Anti-Human IgG(H+L)



www.ptglab.com

Catalog Number:SA00002-7

Information

Catalog Number:	SA00002-7	Reactivity:	Human
Host:	Donkey	Physical State:	Liquid
Applications:	ELISA, WB	Conjugation:	Alkaline phosphatase

Recommended Dilutions

1:1000 – 1:3000 for ELISA and Western blotting

Preservative

0.01% Sodium Azide

Stabilizer

5 mg/ml BSA

Safety Notes

This product is for research use only, not for diagnostic or therapeutic use.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
0.01 M Tris-HCl, 0.25 M NaCl, 50% glycerol, 5 mg/ml BSA, 0.01% Sodium Azide, pH 8.0.
Aliquoting is unnecessary for -20°C storage

Purity

The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.

For technical support and original validation data for this product please contact:

T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data