

For Research Use Only.
Not For Use In Diagnostics.

Fluorescein (FITC)–conjugated Affinipure Goat Anti-Rabbit IgG(H+L)



Catalog Number:SA00003-2

Information

Catalog Number:

SA00003-2

Host:

Goat

Applications:

IF, FC

Reactivity:

Rabbit

Physical State:

Liquid

Conjugation:

FITC

Recommended Dilutions

1:100-1:1000 for most applications

Preservative

0.02% Sodium Azide

Fluorophore

Fluorescein-5-isothiocyanate (FITC) Amax=492nm;Emax=520nm

Fluorophore Protein

9.0 ug/mg; 3.4 moles FITC per mole IgG

Stabilizer

5 mg/ml BSA

Safety Notes

This product is for research use only,not for diagnostic or therapeutic use.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

0.01 M Sodium phosphate, 0.25 M NaCl, 50% glycerol, 5 mg/ml BSA, 0.02% Sodium Azide, pH 7.6.

Aliquoting is unnecessary for -20°C storage

Purity

The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.

For technical support and original validation data for this product please
contact:

T: 4006900926

E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech
Group brand and is not available to purchase from any
other manufacturer.

