

For Research Use Only.
Not For Use In Diagnostics.

Fluorescein (FITC)-conjugated Affinipure Rabbit Anti-Goat IgG(H+L)



www.ptglab.com

Catalog Number: SA00003-4

Information

Catalog Number:	SA00003-4	Reactivity:	Goat
Host:	Rabbit	Physical State:	Liquid
Applications:	IF, FC	Conjugation:	FITC

Recommended Dilutions

1:20-1:100 for most applications

Preservative

0.02% Sodium Azide

Fluorophore

Fluorescein-5-isothiocyanate (FITC) Amax=492nm;Emax=520nm

Fluorophore Protein

9.0 ug/mg; 2.9 moles FITC per mole IgG

Stabilizer

5 mg/ml BSA

Safety Notes

This product is for research use only, not for diagnostic or therapeutic use.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
0.01 M Sodium phosphate, 0.25 M NaCl, 50% glycerol, 5 mg/ml BSA, 0.02% Sodium Azide, pH 7.6.
Aliquoting is unnecessary for -20°C storage

Purity

The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data