

For Research Use Only.  
Not For Use In Diagnostics.

# Biotin-conjugated Affinipure Goat Anti-Rabbit IgG(H+L)

Catalog Number:SA00004-2



## Information

### Catalog Number:

SA00004-2

### Host:

Goat

### Applications:

ELISA, WB, IHC/ICC

### Reactivity:

Rabbit

### Physical State:

Liquid

### Conjugation:

Biotin

## Recommended Dilutions

1:100 – 1:1000 for enzyme immunohistochemistry on tissue sections  
1:4000-1:80000 for ELISA and Western blotting with chromogenic substrates

## Preservative

0.01% Sodium Azide

## Stabilizer

5 mg/ml BSA

## Safety Notes

This product is for research use only, not for diagnostic or therapeutic use.

## Storage

### Storage:

Store at -20°C. Stable for one year after shipment.

### Storage Buffer:

0.01 M Sodium phosphate, 0.25 M NaCl, 50% glycerol, 5 mg/ml BSA, 0.01% Sodium Azide, pH 7.6.

Aliquoting is unnecessary for -20°C storage

## Purity

The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

