

For Research Use Only.
Not For Use In Diagnostics.

Texas Red-conjugated Affinipure Donkey Anti-Goat IgG(H+L)

Catalog Number:SA00005-3



www.ptglab.com

Information

Catalog Number:

SA00005-3

Host:

Donkey

Applications:

IF, FC

Reactivity:

Goat

Physical State:

Liquid

Conjugation:

Texas Red

Recommended Dilutions

1:20-1:100 for most applications

Preservative

0.02% Sodium Azide

Fluorophore

Texas Red sulfonyl chloride (TRSC) $A_{max}=596nm, E_{max}=620nm$

Stabilizer

5 mg/ml BSA

Absorbance Ratio

A596/A280=0.57

Safety Notes

This product is for research use only, not for diagnostic or therapeutic use.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6

Aliquoting is unnecessary for -20°C storage

Purity

The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data