

For Research Use Only.  
Not For Use In Diagnostics.

# Rhodamine (TRITC)-conjugated Donkey Anti-Human IgG(H+L)

Catalog Number:SA00007-11



[www.ptglab.com](http://www.ptglab.com)

## Information

Catalog Number:

SA00007-11

Host:

Donkey

Applications:

IF, FC

Reactivity:

Human

Physical State:

Liquid

Conjugation:

Rhodamine

## Recommended Dilutions

1:20-1:100 for most applications

## Preservative

0.02% Sodium Azide

## Fluorophore

Tetramethyl Rhodamine Isothiocyanate (TRITC)  $A_{max}=550nm, E_{max}=570nm$

## Stabilizer

5 mg/ml BSA

## Absorbance Ratio

A550/A280=0.55

## Safety Notes

This product is for research use only, not for diagnostic or therapeutic use.

## Storage

**Storage:**

Store at  $-20^{\circ}C$ . Stable for one year after shipment.

**Storage Buffer:**

0.01 M Sodium phosphate, 0.25 M NaCl, 50% glycerol, 5 mg/ml BSA, 0.02% Sodium Azide, pH 7.6.

Aliquoting is unnecessary for  $-20^{\circ}C$  storage

## Purity

The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com) W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

---

## Selected Validation Data