

For Research Use Only.
Not For Use In Diagnostics.

R-PE-conjugated AffiniPure F(ab')₂ Fragment Donkey Anti-Mouse IgG(H+L)



www.ptglab.com

Catalog Number: SA00008-9

Information

Catalog Number:

SA00008-9

Host:

Donkey

Applications:

IF, FC

Reactivity:

Mouse

Physical State:

Liquid

Conjugation:

R-PE

Recommended Dilutions

1:20-1:100 for most applications

Preservative

0.02% Sodium Azide

Fluorophore

Phycoerythrin A_{max}=490nm, 545nm and 566nm; E_{max}=580nm

Stabilizer

5 mg/ml BSA

Safety Notes

This product is for research use only, not for diagnostic or therapeutic use.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

0.01 M Sodium phosphate, 0.25 M NaCl, 50% glycerol, 5 mg/ml BSA, 0.02% Sodium Azide, pH 7.6.

Aliquoting is unnecessary for -20°C storage

Purity

The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data