

For Research Use Only.
Not For Use In Diagnostics.

AMCA-conjugated Affinipure Rabbit Anti-Goat IgG(H+L)

Catalog Number:SA00010-4



www.ptglab.com

Information

Catalog Number:

SA00010-4

Host:

Rabbit

Applications:

IF, FC

Reactivity:

Goat

Physical State:

Liquid

Conjugation:

AMCA

Recommended Dilutions

1:20-1:100 for most applications

Preservative

0.02% Sodium Azide

Fluorophore

7-Amino-4-methylcoumarin-3-acetic acid $A_{max}=350nm$, $E_{max}=450nm$

Stabilizer

5 mg/ml BSA

Absorbance Ratio

A350/A280=0.66

Safety Notes

This product is for research use only, not for diagnostic or therapeutic use.

Storage

Storage:

Store at $-20^{\circ}C$. Stable for one year after shipment.

Storage Buffer:

0.01 M Sodium phosphate, 0.25 M NaCl, 50% glycerol, 5 mg/ml BSA, 0.02% Sodium Azide, pH 7.6.

Aliquoting is unnecessary for $-20^{\circ}C$ storage

Purity

The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data