

For Research Use Only.  
Not For Use In Diagnostics.

# HRP-conjugated Affinipure Goat Anti-Mouse IgM, $\mu$ Chain Specific

Catalog Number:SA00012-6



[www.ptglab.com](http://www.ptglab.com)

## Information

**Catalog Number:**

SA00012-6

**Host:**

Goat

**Applications:**

ELISA, WB, Dot blot

**Reactivity:**

Mouse IgM

**Physical State:**

Liquid

**Conjugation:**

HRP

## Recommended Dilutions

1:2000-1:40,000 for western blotting with ECL substrates

1:1000-1:20,000 for ELISA and Western blotting with chromogenic substrates

## Stabilizer

3 mg/ml BSA

## Safety Notes

This product is for research use only, not for diagnostic or therapeutic use.

## Storage

**Storage:**

Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**

0.01 M Sodium phosphate, 0.25 M NaCl, 50% glycerol, 3 mg/ml BSA, pH 7.6.

Aliquoting is unnecessary for -20°C storage

## Purity

The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com) W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data