

For Research Use Only.  
Not For Use In Diagnostics.

# CoraLite594 – conjugated Donkey Anti-Mouse IgG(H+L)

Catalog Number:SA00013-7



[www.ptglab.com](http://www.ptglab.com)

## Information

**Catalog Number:**

SA00013-7

**Host:**

Donkey

**Applications:**

IF, FC

**Reactivity:**

Mouse (Minimal Cross Reactivity:

Bovine, Chicken, Goat, Guinea Pig,

Syrian Hamster, Horse, Human,

Rabbit, Sheep Serum Proteins)

**Physical State:**

Liquid

**Conjugation:**

CoraLite®594

## Recommended Dilutions

1:100 – 1:500

## Preservative

0.02% Sodium Azide

## Fluorophore

CL594 Fluorescent Dye  
Amax=593 nm, Emax=614 nm

## Stabilizer

5 mg/ml BSA

## Safety Notes

This product is for research use only, not for diagnostic or therapeutic use.

## Storage

**Storage:**

Store at -20°C. Stable for one year after shipment.

**Storage Buffer:**

0.01 M Sodium phosphate, 0.25 M NaCl, 50% glycerol, 5 mg/ml BSA, 0.02% Sodium Azide, pH 7.6.

Aliquoting is unnecessary for -20°C storage

## Purity

The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com) W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data