

For Research Use Only.
Not For Use In Diagnostics.

CoraLite594 – conjugated Donkey Anti-Rabbit IgG(H+L)

Catalog Number:SA00013-8



www.ptglab.com

Information

Catalog Number:

SA00013-8

Host:

Donkey

Applications:

IF, FC

Reactivity:

Rabbit (Minimal Cross Reactivity:
Bovine, Chicken, Goat, Guinea Pig,
Syrian Hamster, Horse, Human,
Mouse, Rat, Sheep Serum Proteins)

Physical State:

Liquid

Conjugation:

CoraLite®594

Recommended Dilutions

1:100 – 1:500

Preservative

0.02% Sodium Azide

Fluorophore

CL594 Fluorescent Dye $\lambda_{max}=593\text{ nm}$, $\lambda_{em}=614\text{ nm}$

Stabilizer

5 mg/ml BSA

Safety Notes

This product is for research use only, not for diagnostic or therapeutic use.

Storage

Storage:

Store at -20°C . Stable for one year after shipment.

Storage Buffer:

0.01 M Sodium phosphate, 0.25 M NaCl, 50% glycerol, 5 mg/ml BSA, 0.02% Sodium Azide, pH 7.6.

Aliquoting is unnecessary for -20°C storage

Purity

The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data